

## Product datasheet for **KN209139RB**

### SOX4 Human Gene Knockout Kit (CRISPR)

#### Product data:

**Product Type:** Knockout Kits (CRISPR)

**Format:** 2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control

**Donor DNA:** RFP-BSD

**Symbol:** SOX4

**Locus ID:** 6659

**Components:** **KN209139G1**, SOX4 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002)  
**KN209139G2**, SOX4 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002)  
**KN209139RBD**, donor DNA containing left and right homologous arms and RFP-BSD functional cassette.  
**GE100003**, scramble sequence in pCas-Guide vector

**Disclaimer:** These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.

**RefSeq:** [NM\\_003107](#)

**UniProt ID:** [Q06945](#)

**Synonyms:** EVI16

**Summary:** This intronless gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins, such as syndecan binding protein (syntenin). The protein may function in the apoptosis pathway leading to cell death as well as to tumorigenesis and may mediate downstream effects of parathyroid hormone (PTH) and PTH-related protein (PTHrP) in bone development. The solution structure has been resolved for the HMG-box of a similar mouse protein. [provided by RefSeq, Jul 2008]



[View online »](#)

## Product images:

