

Product datasheet for **KN209243BN**

ZNF281 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	ZNF281
Locus ID:	23528
Components:	KN209243G1 , ZNF281 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN209243G2 , ZNF281 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN209243BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001281293 , NM_001281294 , NM_012482
UniProt ID:	Q9Y2X9
Synonyms:	ZBP-99; ZNP-99
Summary:	Transcription repressor that plays a role in regulation of embryonic stem cells (ESCs) differentiation. Required for ESCs differentiation and acts by mediating autorepression of NANOG in ESCs: binds to the NANOG promoter and promotes association of NANOG protein to its own promoter and recruits the NuRD complex, which deacetylates histones. Not required for establishment and maintenance of ESCs (By similarity). Represses the transcription of a number of genes including GAST, ODC1 and VIM. Binds to the G-rich box in the enhancer region of these genes.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

