

Product datasheet for **KN209340LP**

SCP3 (SYCP3) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	SCP3
Locus ID:	50511
Components:	KN209340G1 , SCP3 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN209340G2 , SCP3 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN209340LPD , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001177948 , NM_001177949 , NM_153694
UniProt ID:	Q8IZU3
Synonyms:	COR1; SCP3; SPGF4
Summary:	This gene encodes an essential structural component of the synaptonemal complex. This complex is involved in synapsis, recombination and segregation of meiotic chromosomes. Mutations in this gene are associated with azoospermia in males and susceptibility to pregnancy loss in females. Alternate splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, May 2010]



[View online »](#)

Product images:

