

Product datasheet for **KN209597BN**

IMP3 (IGF2BP3) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	IMP3
Locus ID:	10643
Components:	KN209597G1 , IMP3 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN209597G2 , IMP3 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN209597BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_006547
UniProt ID:	O00425
Synonyms:	CT98; IMP-3; IMP3; KOC; KOC1; VICKZ3
Summary:	The protein encoded by this gene is primarily found in the nucleolus, where it can bind to the 5' UTR of the insulin-like growth factor II leader 3 mRNA and may repress translation of insulin-like growth factor II during late development. The encoded protein contains several KH domains, which are important in RNA binding and are known to be involved in RNA synthesis and metabolism. A pseudogene exists on chromosome 7, and there are putative pseudogenes on other chromosomes. [provided by RefSeq, Jul 2008]



[View online »](#)

Product images:
