

## Product datasheet for **KN209832BN**

### Aryl hydrocarbon Receptor (AHR) Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	Aryl hydrocarbon Receptor
Locus ID:	196
Components:	<b>KN209832G1</b> , Aryl hydrocarbon Receptor gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN209832G2</b> , Aryl hydrocarbon Receptor gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN209832BND</b> , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_001621</a>
UniProt ID:	<a href="#">P35869</a>
Synonyms:	bHLHe76
Summary:	The protein encoded by this gene is a ligand-activated helix-loop-helix transcription factor involved in the regulation of biological responses to planar aromatic hydrocarbons. This receptor has been shown to regulate xenobiotic-metabolizing enzymes such as cytochrome P450. Before ligand binding, the encoded protein is sequestered in the cytoplasm; upon ligand binding, this protein moves to the nucleus and stimulates transcription of target genes. [provided by RefSeq, Sep 2015]



[View online »](#)

Product images:

