

Product datasheet for **KN209855LP**

SLC25A24 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	SLC25A24
Locus ID:	29957
Components:	KN209855G1 , SLC25A24 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN209855G2 , SLC25A24 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN209855LPD , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_013386 , NM_213651
UniProt ID:	Q6NUK1
Synonyms:	APC1; SCAMC-1
Summary:	This gene encodes a carrier protein that transports ATP-Mg exchanging it for phosphate. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2012]



[View online »](#)

Product images:

