

## Product datasheet for **KN210171BN**

### Calcitonin receptor (CALCR) Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	Calcitonin receptor
Locus ID:	799
Components:	<b>KN210171G1</b> , Calcitonin receptor gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN210171G2</b> , Calcitonin receptor gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN210171BND</b> , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_001164737</a> , <a href="#">NM_001164738</a> , <a href="#">NM_001742</a>
UniProt ID:	<a href="#">P30988</a>
Synonyms:	CRT; CT-R; CTR; CTR1
Summary:	This gene encodes a high affinity receptor for the peptide hormone calcitonin and belongs to a subfamily of seven transmembrane-spanning G protein-coupled receptors. The encoded protein is involved in maintaining calcium homeostasis and in regulating osteoclast-mediated bone resorption. Polymorphisms in this gene have been associated with variations in bone mineral density and onset of osteoporosis. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2009]



[View online »](#)

## Product images:

