

Product datasheet for **KN210211RB**

12 Lipoxygenase (ALOX12) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	12 Lipoxygenase
Locus ID:	239
Components:	KN210211G1 , 12 Lipoxygenase gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN210211G2 , 12 Lipoxygenase gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN210211RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_000697
UniProt ID:	P18054
Synonyms:	12-LOX; 12S-LOX; LOG12
Summary:	This gene encodes a member of the lipoxygenase family of proteins. The encoded enzyme acts on different polyunsaturated fatty acid substrates to generate bioactive lipid mediators including eicosanoids and lipoxins. The encoded enzyme and its reaction products have been shown to regulate platelet function. Elevated expression of this gene has been observed in pancreatic islets derived from human diabetes patients. Allelic variants in this gene may be associated with susceptibility to toxoplasmosis. Multiple pseudogenes of this gene have been identified in the human genome. [provided by RefSeq, Aug 2017]



[View online »](#)

Product images:

