

## Product datasheet for **KN211218LP**

### Tyrosine Hydroxylase (TH) Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	Tyrosine Hydroxylase
Locus ID:	7054
Components:	<b>KN211218G1</b> , Tyrosine Hydroxylase gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN211218G2</b> , Tyrosine Hydroxylase gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN211218LPD</b> , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_000360</a> , <a href="#">NM_199292</a> , <a href="#">NM_199293</a>
UniProt ID:	<a href="#">P07101</a>
Synonyms:	DYT5b; DYT14; TYH
Summary:	The protein encoded by this gene is involved in the conversion of tyrosine to dopamine. It is the rate-limiting enzyme in the synthesis of catecholamines, hence plays a key role in the physiology of adrenergic neurons. Mutations in this gene have been associated with autosomal recessive Segawa syndrome. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Jul 2008]



[View online »](#)

## Product images:

