

Product datasheet for **KN211582RB**

BTK Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	BTK
Locus ID:	695
Components:	KN211582G1 , BTK gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN211582G2 , BTK gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN211582RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_000061 , NM_001287344 , NM_001287345
UniProt ID:	Q06187
Synonyms:	AGMX1; AT; ATK; BPK; IMD1; PSCTK1; XLA
Summary:	The protein encoded by this gene plays a crucial role in B-cell development. Mutations in this gene cause X-linked agammaglobulinemia type 1, which is an immunodeficiency characterized by the failure to produce mature B lymphocytes, and associated with a failure of Ig heavy chain rearrangement. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2013]



[View online »](#)

Product images:

