

## Product datasheet for **KN211665RB**

### KSR2 Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	KSR2
Locus ID:	283455
Components:	<b>KN211665G1</b> , KSR2 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN211665G2</b> , KSR2 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN211665RBD</b> , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_173598</a>
UniProt ID:	<a href="#">Q6VAB6</a>
Synonyms:	FLJ25965
Summary:	Location-regulated scaffold connecting MEK to RAF. Has very low protein kinase activity and can phosphorylate MAP2K1 at several Ser and Thr residues with very low efficiency (in vitro). Interaction with BRAF enhances KSR2-mediated phosphorylation of MAP2K1 (in vitro). Blocks MAP3K8 kinase activity and MAP3K8-mediated signaling. Acts as a negative regulator of MAP3K3-mediated activation of ERK, JNK and NF-kappa-B pathways, inhibiting MAP3K3-mediated interleukin-8 production.[UniProtKB/Swiss-Prot Function]



[View online »](#)

## Product images:

