

Product datasheet for **KN213508RB**

PNMA8B Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	PNMA8B
Locus ID:	57469
Components:	KN213508G1 , PNMA8B gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN213508G2 , PNMA8B gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN213508RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.

RefSeq:	NM_020709
UniProt ID:	Q9ULN7
Synonyms:	FLJ42730; KIAA1183



[View online »](#)

Product images:

