

## Product datasheet for **KN215023BN**

### IRE1 (ERN1) Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	IRE1
Locus ID:	2081
Components:	<p><b>KN215023G1</b>, IRE1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002)</p> <p><b>KN215023G2</b>, IRE1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002)</p> <p><b>KN215023BND</b>, donor DNA containing left and right homologous arms and mBFP-Neo functional cassette.</p> <p><b>GE100003</b>, scramble sequence in pCas-Guide vector</p>
Disclaimer:	<p>These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.</p>
RefSeq:	<u><a href="#">NM_001433</a></u> , <u><a href="#">NM_152461</a></u>
UniProt ID:	<u><a href="#">O75460</a></u>
Synonyms:	hIRE1p; IRE1; IRE1a; IRE1P
Summary:	<p>This gene encodes the transmembrane protein kinase inositol-requiring enzyme 1. The encoded protein contains two functional catalytic domains, a serine/threonine-protein kinase domain and an endoribonuclease domain. This protein functions as a sensor of unfolded proteins in the endoplasmic reticulum (ER) and triggers an intracellular signaling pathway termed the unfolded protein response (UPR). The UPR is an ER stress response that is conserved from yeast to mammals and activates genes involved in degrading misfolded proteins, regulating protein synthesis and activating molecular chaperones. This protein specifically mediates the splicing and activation of the stress response transcription factor X-box binding protein 1. [provided by RefSeq, Aug 2017]</p>



[View online »](#)

## Product images:

