

Product datasheet for **KN215065BN**

NLRP7 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	NLRP7
Locus ID:	199713
Components:	KN215065G1 , NLRP7 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN215065G2 , NLRP7 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN215065BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001127255 , NM_139176 , NM_206828
UniProt ID:	Q8WX94
Synonyms:	CLR19.4; HYDM; NALP7; NOD12; PAN7; PYPAF3
Summary:	This gene encodes a member of the NACHT, leucine rich repeat, and PYD containing (NLRP) protein family. It has an N-terminal pyrin domain, followed by a NACHT domain, a NACHT-associated domain (NAD), and a C-terminal leucine-rich repeat (LRR) region. NLRP proteins are implicated in the activation of proinflammatory caspases through multiprotein complexes called inflammasomes. This gene may act as a feedback regulator of caspase-1-dependent interleukin 1-beta secretion. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]



[View online »](#)

Product images:

