

## Product datasheet for **KN216942LP**

### Neogenin (NEO1) Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	Neogenin
Locus ID:	4756
Components:	<b>KN216942G1</b> , Neogenin gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN216942G2</b> , Neogenin gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN216942LPD</b> , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_001172623</a> , <a href="#">NM_001172624</a> , <a href="#">NM_002499</a>
UniProt ID:	<a href="#">Q92859</a>
Synonyms:	IGDCC2; NGN; NTN1R2
Summary:	This gene encodes a cell surface protein that is a member of the immunoglobulin superfamily. The encoded protein consists of four N-terminal immunoglobulin-like domains, six fibronectin type III domains, a transmembrane domain and a C-terminal internal domain that shares homology with the tumor suppressor candidate gene DCC. This protein may be involved in cell growth and differentiation and in cell-cell adhesion. Defects in this gene are associated with cell proliferation in certain cancers. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Feb 2010]



[View online »](#)

## Product images:

