

## Product datasheet for **KN217426RB**

### FGF2 Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	FGF2
Locus ID:	2247
Components:	<p><b>KN217426G1</b>, FGF2 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002)</p> <p><b>KN217426G2</b>, FGF2 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002)</p> <p><b>KN217426RBD</b>, donor DNA containing left and right homologous arms and RFP-BSD functional cassette.</p> <p><b>GE100003</b>, scramble sequence in pCas-Guide vector</p>
Disclaimer:	<p>These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.</p>
RefSeq:	<u><a href="#">NM_002006</a></u> , <u><a href="#">NM_001361665</a></u>
UniProt ID:	<u><a href="#">P09038</a></u>
Synonyms:	BFGF; FGF-2; FGFB; HBGF-2
Summary:	<p>The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members bind heparin and possess broad mitogenic and angiogenic activities. This protein has been implicated in diverse biological processes, such as limb and nervous system development, wound healing, and tumor growth. The mRNA for this gene contains multiple polyadenylation sites, and is alternatively translated from non-AUG (CUG) and AUG initiation codons, resulting in five different isoforms with distinct properties. The CUG-initiated isoforms are localized in the nucleus and are responsible for the intracrine effect, whereas, the AUG-initiated form is mostly cytosolic and is responsible for the paracrine and autocrine effects of this FGF. [provided by RefSeq, Jul 2008]</p>



[View online »](#)

## Product images:

