

## Product datasheet for **KN218078BN**

### Argonaute 2 (AGO2) Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	Argonaute 2
Locus ID:	27161
Components:	<b>KN218078G1</b> , Argonaute 2 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN218078G2</b> , Argonaute 2 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN218078BND</b> , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_001164623</a> , <a href="#">NM_012154</a>
UniProt ID:	<a href="#">Q9UKV8</a>
Synonyms:	EIF2C2; Q10
Summary:	This gene encodes a member of the Argonaute family of proteins which play a role in RNA interference. The encoded protein is highly basic, and contains a PAZ domain and a PIWI domain. It may interact with dicer1 and play a role in short-interfering-RNA-mediated gene silencing. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]



[View online »](#)

Product images:

