

Product datasheet for KN218078RB

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Argonaute 2 (AGO2) Human Gene Knockout Kit (CRISPR)

Product data:

Product Type: Knockout Kits (CRISPR)

Format: 2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control

Donor DNA: RFP-BSD
Symbol: Argonaute 2

Locus ID: 27161

Components: KN218078G1, Argonaute 2 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002)

KN218078G2, Argonaute 2 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) **KN218078RBD**, donor DNA containing left and right homologous arms and RFP-BSD

functional cassette.

GE100003, scramble sequence in pCas-Guide vector

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS. The kit is

designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the

experimental process.

RefSeq: <u>NM 001164623</u>, <u>NM 012154</u>

UniProt ID: Q9UKV8

Synonyms: EIF2C2; Q10

Summary: This gene encodes a member of the Argonaute family of proteins which play a role in RNA

interference. The encoded protein is highly basic, and contains a PAZ domain and a PIWI domain. It may interact with dicer1 and play a role in short-interfering-RNA-mediated gene silencing. Multiple transcript variants encoding different isoforms have been found for this


gene. [provided by RefSeq, Sep 2009]





Product images:

Donor Vector Edited Chromosome



RFP, Luc, and mBFP will be under native gene promoter