

## Product datasheet for **KN218134LP**

### SIRT1 Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	SIRT1
Locus ID:	23411
Components:	<b>KN218134G1</b> , SIRT1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN218134G2</b> , SIRT1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN218134LPD</b> , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_001142498</a> , <a href="#">NM_001314049</a> , <a href="#">NM_012238</a>
UniProt ID:	<a href="#">Q96EB6</a>
Synonyms:	SIR2; SIR2alpha; SIR2L1
Summary:	This gene encodes a member of the sirtuin family of proteins, homologs to the yeast Sir2 protein. Members of the sirtuin family are characterized by a sirtuin core domain and grouped into four classes. The functions of human sirtuins have not yet been determined; however, yeast sirtuin proteins are known to regulate epigenetic gene silencing and suppress recombination of rDNA. Studies suggest that the human sirtuins may function as intracellular regulatory proteins with mono-ADP-ribosyltransferase activity. The protein encoded by this gene is included in class I of the sirtuin family. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2008]



[View online »](#)

Product images:

