

Product datasheet for **KN220046LP**

RAD54 (RAD54L) Human Gene Knockout Kit (CRISPR)

Product data:

| | |
|---------------|---|
| Product Type: | Knockout Kits (CRISPR) |
| Format: | 2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control |
| Donor DNA: | Luciferase-Puro |
| Symbol: | RAD54 |
| Locus ID: | 8438 |
| Components: | KN220046G1 , RAD54 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN220046G2 , RAD54 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN220046LPD , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. GE100003 , scramble sequence in pCas-Guide vector |
| Disclaimer: | These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process. |
| RefSeq: | NM_001142548 , NM_003579 , NM_001370766 |
| UniProt ID: | Q92698 |
| Synonyms: | hHR54; HR54; hRAD54; RAD54A |
| Summary: | The protein encoded by this gene belongs to the DEAD-like helicase superfamily, and shares similarity with <i>Saccharomyces cerevisiae</i> Rad54, a protein known to be involved in the homologous recombination and repair of DNA. This protein has been shown to play a role in homologous recombination related repair of DNA double-strand breaks. The binding of this protein to double-strand DNA induces a DNA topological change, which is thought to facilitate homologous DNA pairing, and stimulate DNA recombination. Alternative splicing results in multiple transcript variants encoding the same protein.[provided by RefSeq, Dec 2008] |



[View online »](#)

Product images:

