

## Product datasheet for **KN221327BN**

### **Snf1lk2 (SIK2) Human Gene Knockout Kit (CRISPR)**

#### **Product data:**

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	Snf1lk2
Locus ID:	23235
Components:	<b>KN221327G1</b> , Snf1lk2 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN221327G2</b> , Snf1lk2 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN221327BND</b> , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_015191</a>
UniProt ID:	<a href="#">Q9H0K1</a>
Synonyms:	LOH11CR11; QIK; SIK-2; SNF1LK2
Summary:	Phosphorylates 'Ser-794' of IRS1 in insulin-stimulated adipocytes, potentially modulating the efficiency of insulin signal transduction. Inhibits CREB activity by phosphorylating and repressing TORCs, the CREB-specific coactivators.[UniProtKB/Swiss-Prot Function]



[View online »](#)

## Product images:

