

## Product datasheet for **KN222150LP**

### CCHCR1 Human Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	CCHCR1
Locus ID:	54535
Components:	<b>KN222150G1</b> , CCHCR1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN222150G2</b> , CCHCR1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN222150LPD</b> , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_001105563</a> , <a href="#">NM_001105564</a> , <a href="#">NM_019052</a>
UniProt ID:	<a href="#">Q8TD31</a>
Synonyms:	C6orf18; HCR; SBP
Summary:	This gene encodes a protein with five coiled-coil alpha-helical rod domains that is thought to act as a regulator of mRNA metabolism through its interaction with mRNA-decapping protein 4. It localizes to P-bodies, the site of mRNA metabolism, with an N-terminus that is required for this subcellular localization, suggesting it is a P-body component. Naturally occurring mutations in this gene are associated with psoriasis. [provided by RefSeq, May 2017]



[View online »](#)

Product images:

