

Product datasheet for **KN222399BN**

DLL1 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	DLL1
Locus ID:	28514
Components:	KN222399G1 , DLL1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN222399G2 , DLL1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN222399BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_005618
UniProt ID:	O00548
Synonyms:	Delta; DELTA1; DL1
Summary:	DLL1 is a human homolog of the Notch Delta ligand and is a member of the delta/serrate/jagged family. It plays a role in mediating cell fate decisions during hematopoiesis. It may play a role in cell-to-cell communication. [provided by RefSeq, Jul 2008]



[View online »](#)

Product images:

