

Product datasheet for **KN223140BN**

LMTK3 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	LMTK3
Locus ID:	114783
Components:	KN223140G1 , LMTK3 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN223140G2 , LMTK3 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN223140BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001080434
Synonyms:	LMR3; PPP1R101; TYKLM3



[View online »](#)

Product images:

