

Product datasheet for **KN224753LP**

B4GALNT4 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	B4GALNT4
Locus ID:	338707
Components:	KN224753G1 , B4GALNT4 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN224753G2 , B4GALNT4 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN224753LPD , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_178537
UniProt ID:	Q76KP1
Synonyms:	FLJ25045
Summary:	Transfers N-acetylgalactosamine (GalNAc) from UDP-GalNAc to N-acetylglucosamine-beta-benzyl with a beta-1,4-linkage to form N,N'-diacetyllactosediamine, GalNAc-beta-1,4-GlcNAc structures in N-linked glycans and probably O-linked glycans.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

