

## Product datasheet for **KN224804BN**

### **GAS1 Human Gene Knockout Kit (CRISPR)**

#### **Product data:**

|               |   |
|---------------|---|
| Product Type: | Knockout Kits (CRISPR)  |
| Format:       | 2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control  |
| Donor DNA:    | mBFP-Neo  |
| Symbol:       | GAS1  |
| Locus ID:     | 2619  |
| Components:   | <b>KN224804G1</b> , GAS1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002)<br><b>KN224804G2</b> , GAS1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002)<br><b>KN224804BND</b> , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette.<br><b>GE100003</b> , scramble sequence in pCas-Guide vector  |
| Disclaimer:   | These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process. |
| RefSeq:       | <a href="#">NM_002048</a>   |
| UniProt ID:   | <a href="#">P54826</a>  |
| Synonyms:     | growth arrest-specific 1; Growth arrest-specific gene-1; OTTHUMP00000021599   |
| Summary:      | Growth arrest-specific 1 plays a role in growth suppression. GAS1 blocks entry to S phase and prevents cycling of normal and transformed cells. Gas1 is a putative tumor suppressor gene. [provided by RefSeq, Jul 2008]  |



[View online »](#)

**Product images:**
