

Product datasheet for **KN233808LP**

PLK5 Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	PLK5
Locus ID:	126520
Components:	KN233808G1 , PLK5 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN233808G2 , PLK5 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN233808LPD , donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001011716 , NM_001243079 , NR_026557
UniProt ID:	Q496M5
Synonyms:	PLK-5; PLK5P; SgK384ps
Summary:	Inactive serine/threonine-protein kinase that plays a role in cell cycle progression and neuronal differentiation.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

