

Product datasheet for **KN301242RB**

Angptl4 Mouse Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	Angptl4
Locus ID:	57875
Components:	KN301242G1 , Angptl4 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN301242G2 , Angptl4 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN301242RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.

RefSeq: [NM_020581](#)

UniProt ID: [Q9Z1P8](#)

Synonyms: Arp4; Bk89; Fiaf; Hfarp; Ng27; Pgar; Pgarg; Pp1158

Summary: Mediates inactivation of the lipoprotein lipase LPL, and thereby plays a role in the regulation of triglyceride clearance from the blood serum and in lipid metabolism (PubMed:15837923, PubMed:17609370, PubMed:29899519). May also play a role in regulating glucose homeostasis and insulin sensitivity (PubMed:15837923, PubMed:29899519). Inhibits proliferation, migration, and tubule formation of endothelial cells and reduces vascular leakage (PubMed:14583458, PubMed:17130448, PubMed:21832056). Upon heterologous expression, inhibits the adhesion of endothelial cell to the extracellular matrix (ECM), and inhibits the reorganization of the actin cytoskeleton, formation of actin stress fibers and focal adhesions in endothelial cells that have adhered to ANGPTL4-containing ECM (in vitro) (By similarity). Depending on context, may modulate tumor-related angiogenesis (Probable). [UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

