

## Product datasheet for **KN304621LP**

### Dll3 Mouse Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 Luciferase-Puro donor, 1 scramble control
Donor DNA:	Luciferase-Puro
Symbol:	Dll3
Locus ID:	13389
Components:	<p><b>KN304621G1</b>, Dll3 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002)</p> <p><b>KN304621G2</b>, Dll3 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002)</p> <p><b>KN304621LPD</b>, donor DNA containing left and right homologous arms and Luciferase-Puro functional cassette.</p> <p><b>GE100003</b>, scramble sequence in pCas-Guide vector</p>
Disclaimer:	<p>These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.</p>
RefSeq:	<a href="#">NM_007866</a>
UniProt ID:	<a href="#">O88516</a>
Synonyms:	pu; pudgy
Summary:	Inhibits primary neurogenesis. May be required to divert neurons along a specific differentiation pathway. Plays a role in the formation of somite boundaries during segmentation of the paraxial mesoderm.[UniProtKB/Swiss-Prot Function]



[View online »](#)

## Product images:

