

Product datasheet for **KN304670BN**

Dnajb3 Mouse Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	Dnajb3
Locus ID:	15504
Components:	KN304670G1 , Dnajb3 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN304670G2 , Dnajb3 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN304670BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_008299
UniProt ID:	O35723
Synonyms:	Hsj3; MSJ-1; Msj1
Summary:	May operate as a co-chaperone of the male germ cell- and haploid stage-specific Hsp70 proteins.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

