

Product datasheet for **KN304882RB**

Dusp6 Mouse Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	Dusp6
Locus ID:	67603
Components:	KN304882G1 , Dusp6 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN304882G2 , Dusp6 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN304882RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_026268
UniProt ID:	Q9DBB1
Synonyms:	1300019I03Rik; MKP-3; MKP3; PYST1
Summary:	Inactivates MAP kinases. Has a specificity for the ERK family (By similarity). Plays an important role in alleviating acute postoperative pain (PubMed:24155322, PubMed:28405172). Necessary for the normal dephosphorylation of the long-lasting phosphorylated forms of spinal MAPK1/3 and MAP kinase p38 induced by peripheral surgery, which drives the resolution of acute postoperative allodynia (PubMed:24155322). Also important for dephosphorylation of MAPK1/3 in local wound tissue, which further contributes to resolution of acute pain (PubMed:28405172).[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

