

## Product datasheet for **KN308871RB**

### **Klhl32 Mouse Gene Knockout Kit (CRISPR)**

#### **Product data:**

<b>Product Type:</b>	Knockout Kits (CRISPR)
<b>Format:</b>	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
<b>Donor DNA:</b>	RFP-BSD
<b>Symbol:</b>	Klhl32
<b>Locus ID:</b>	212390
<b>Components:</b>	<b>KN308871G1</b> , Klhl32 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN308871G2</b> , Klhl32 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN308871RBD</b> , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
<b>Disclaimer:</b>	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
<b>RefSeq:</b>	<a href="#">NM_001033531</a> , <a href="#">NM_001163020</a> , <a href="#">NM_001355618</a>
<b>Synonyms:</b>	6430524H05Rik; D4Ertd389e; Gm1356; mKIAA1900



[View online »](#)

## Product images:

