

## Product datasheet for **KN309105RB**

### Lair1 Mouse Gene Knockout Kit (CRISPR)

#### Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	Lair1
Locus ID:	52855
Components:	<b>KN309105G1</b> , Lair1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) <b>KN309105G2</b> , Lair1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) <b>KN309105RBD</b> , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. <b>GE100003</b> , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	<a href="#">NM_001113474</a> , <a href="#">NM_001302675</a> , <a href="#">NM_001302676</a> , <a href="#">NM_001302677</a> , <a href="#">NM_001302681</a> , <a href="#">NM_001302682</a> , <a href="#">NM_001302683</a> , <a href="#">NM_178611</a>
UniProt ID:	<a href="#">Q8BG84</a>
Synonyms:	5133400O11Rik; BB115266; D7Bwg0421e; Lair-1



[View online »](#)

**Summary:**

Functions as an inhibitory receptor that plays a constitutive negative regulatory role on cytolytic function of natural killer (NK) cells, B-cells and T-cells. Activation by Tyr phosphorylation results in recruitment and activation of the phosphatases PTPN6 and PTPN11. It also reduces the increase of intracellular calcium evoked by B-cell receptor ligation. May also play its inhibitory role independently of SH2-containing phosphatases. Modulates cytokine production in CD4+ T-cells, down-regulating IL2 and IFNG production while inducing secretion of transforming growth factor beta. Down-regulates also IgG and IgE production in B-cells as well as IL8, IL10 and TNF secretion. Inhibits proliferation and induces apoptosis in myeloid leukemia cell lines as well as prevents nuclear translocation of NF-kappa-B p65 subunit/RELA and phosphorylation of I-kappa-B alpha/CHUK in these cells. Inhibits the differentiation of peripheral blood precursors towards dendritic cells (By similarity).[UniProtKB/Swiss-Prot Function]

**Product images:**
