

Product datasheet for **KN309965RB**

Meis1 Mouse Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	Meis1
Locus ID:	17268
Components:	KN309965G1 , Meis1 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN309965G2 , Meis1 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN309965RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001193271 , NM_010789
UniProt ID:	Q60954
Synonyms:	C530044H18Rik; Evi8
Summary:	Acts as a transcriptional regulator of PAX6. Also acts as a transcriptional activator of PF4 in complex with PBX1 or PBX2. Required for hematopoiesis, megakaryocyte lineage development and vascular patterning. May function as a cofactor for HOXA7 and HOXA9 in the induction of myeloid leukemias.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

