

Product datasheet for **KN313408RB**

Plaur Mouse Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	Plaur
Locus ID:	18793
Components:	KN313408G1 , Plaur gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN313408G2 , Plaur gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN313408RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_011113
UniProt ID:	P35456
Synonyms:	Cd87; u-PAR; uPAR
Summary:	Acts as a receptor for urokinase plasminogen activator. Plays a role in localizing and promoting plasmin formation. Mediates the proteolysis-independent signal transduction activation effects of U-PA.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

