

Product datasheet for **KN316145BN**

Slc4a5 Mouse Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	Slc4a5
Locus ID:	232156
Components:	KN316145G1 , Slc4a5 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN316145G2 , Slc4a5 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN316145BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001166067
Synonyms:	C330006M16; C330016K18Rik



[View online »](#)

Product images:

