

Product datasheet for **KN317729BN**

Tmem143 Mouse Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	Tmem143
Locus ID:	70209
Components:	KN317729G1 , Tmem143 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN317729G2 , Tmem143 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN317729BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_144801
UniProt ID:	Q8VD26
Synonyms:	2310076O21Rik; AI481604



[View online »](#)

Product images:

