

Product datasheet for **KN318125BN**

Traf3ip2 Mouse Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 mBFP-Neo donor, 1 scramble control
Donor DNA:	mBFP-Neo
Symbol:	Traf3ip2
Locus ID:	103213
Components:	KN318125G1 , Traf3ip2 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN318125G2 , Traf3ip2 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN318125BND , donor DNA containing left and right homologous arms and mBFP-Neo functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_134000
UniProt ID:	Q8N7N6
Synonyms:	Act1; AI429613; CIKS
Summary:	E3 ubiquitin ligase that catalyzes 'Lys63'-linked polyubiquitination of target protein, enhancing protein-protein interaction and cell signaling (By similarity). Transfers ubiquitin from E2 ubiquitin-conjugating enzyme UBE2V1-UBE2N to substrate protein (By similarity). Essential adapter molecule in IL17A-mediated signaling (PubMed:19825828). Upon IL17A stimulation, interacts with IL17RA and IL17RC receptor chains through SEFIR domains and catalyzes 'Lys63'-linked polyubiquitination of TRAF6, leading to TRAF6-mediated activation of NF-kappa-B and MAPkinase pathways (PubMed:19825828).[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

