

Product datasheet for **KN319279RB**

Vsig8 Mouse Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 RFP-BSD donor, 1 scramble control
Donor DNA:	RFP-BSD
Symbol:	Vsig8
Locus ID:	240916
Components:	KN319279G1 , Vsig8 gRNA vector 1 in pCas-Guide CRISPR vector (GE100002) KN319279G2 , Vsig8 gRNA vector 2 in pCas-Guide CRISPR vector (GE100002) KN319279RBD , donor DNA containing left and right homologous arms and RFP-BSD functional cassette. GE100003 , scramble sequence in pCas-Guide vector
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_177723 , NR_027644
UniProt ID:	Q6P3A4
Synonyms:	A030011M19; EG240916



[View online »](#)

Product images:

