

Product datasheet for **KN419883**

SPECCIL Human Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 linear donor
Donor DNA:	EF1a-GFP-P2A-Puro
Symbol:	SPECCIL
Locus ID:	23384



Components:

KN419883G1, SPECCIL gRNA vector 1 in pCas-Guide CRISPR vector (GE100002)

KN419883G2, SPECCIL gRNA vector 2 in pCas-Guide CRISPR vector (GE100002)

KN419883D, Linear donor DNA containing LoxP-EF1A-tGFP-P2A-Puro-LoxP:

The sequence below is cassette sequence only. The linear donor DNA also contains proprietary target sequence.

LoxP-EF1A-tGFP-P2A-Puro-LoxP (2739 bp)

ATAACTTCGT ATAATGTATG CTATACGAAG TTATCGTGAG GCTCCGGTGC CCGTCAGTGG GCAGAGCGCA CATCGCCAC
AGTCCCGCAG AAGTTGGGGG GAGGGGTCCG CAATTGAACC GGTGCCTAGA GAAGGTGGCG CGGGGTAAC TGGGAAAGTG
ATGTCGTGTA CTGGCTCCGC CTTTTTCCCG AGGGTGGGG AGAACCGTAT ATAAGTCAG TAGTCGCCGT GAACGTTCTT
TTTCGCAACG GGTTCGCCG CAGAACACAG GTAAGTGCCG TGTGTGGTTC CCGCGGGCCT GGCCTCTTTA CGGGTTATGG
CCCTTGCGTG CCTTGAATTA CTTCCACCTG GCTGCAGTAC GTGATTCTTG ATCCCGAGCT TCGGGTTGGA AGTGGGTGGG
AGAGTTCGAG GCCTTGCGCT TAAGGAGCCC CTTCGCCTCG TGCTTGAGTT GAGGCTGGC CTGGGCGCTG GGGCCGCCG
GTGCGAATCT GGTGGCACCT TCGCGCTGT CTCGCTGCT TCGATAAGTC TCTAGCCATT TAAATTTTTT GATGACCTGC
TGCGACGCTT TTTTCTGGC AAGATAGTCT TGTAATGCG GGCCAAGATC TGCACACTGG TATTTCGGTT TTTGGGGCCG
CGGGCGGCGA CGGGGCCCGT GCGTCCCAGC GCACATGTTT GGCGAGGCGG GGCTGCGAG CGCGGCCACC GAGAATCGGA
CGGGGTAGT CTCAAGCTGG CCGGCTGCT CTGGTGCCTG GCCTCGCGCC GCCGTGTATC GCGCCGCCCT GGGCGGCAAG
GCTGGCCCGG TCGGCACCAG TTGCGTGAGC GGAAGATGG CCGCTCCCG GCCCTGCTGC AGGGAGCTCA AAATGGAGGA
CGCGGCGCTC GGGAGAGCGG GCGGGTGAAG CACCCACACA AAGGAAAAGG GCCTTTCCTG CCTCAGCCGT CGCTTCATGT
GACTCCACGG AGTACCGGGC GCCGTCCAGG CACCTCGATT AGTTCCTGAG CTTTGGAGT ACCTCGTCTT TAGGTTGGGG
GGAGGGGTTT TATGCGATGG AGTTTCCCA CACTGAGTGG GTGGAGACTG AAGTTAGGCC AGCTTGGCAC TTGATGTAAT
TCTCCTGGA ATTTGCCCTT TTTGAGTTG GATCTTGTT CATTCTCAAG CCTCAGACAG TGTTTCAAAG TTTTCTCTT
CCATTTACAG TGTCTGAAT GGAGAGCGAC GAGAGCGGCG TGCCCGCCAT GGAGATCGAG TGCCGCATCA CCGGCACCTC
GAACGCGCTG GAGTTCGAGC TGGTGGGCGG CGGAGAGGCG ACCCCGAGC AGGGCCGCAT GACCAACAAG ATGAAGAGCA
CCAAAGGCGC CCTGACCTTC AGCCCTACC TGCTGAGCCA CGTGATGGG TACGGCTTCT ACCACTTCGG CACCTACCCC
AGCGGCTACG AGAACCCCTT CCTGCACGCC ATCAACAACG GCGGCTACAC CAACACCCG ATCGAGAAGT ACGAGGACGG
CGGCGTGCTG CACGTGAGCT TCAGCTACCG CTACGAGGCC GGCCGCGTGA TCGCGCACTT CAAGGTGATG GGCACCGGCT
TCCCCGAGGA CAGCGTGATC TTCACCGACA AGATCATCCG CAGCAACGCC ACCGTGGAGC ACCTGCACCC CATGGGCGAT
AACGATCTGG ATGGCAGCTT CACCCGCACC TTCAGCCTGC GCGACGGCGG CTAACAGC TCCGTGGTGG ACAGCCACAT
GCATTTCAAG AGCGCCATCC ACCCCAGCAT CCTGCAGAAC GGGGGCCCCA TGTTGCCTT CCGCCGCGTG GAGGAGGATC
ACAGCAACAC CGAGCTGGG ATCGTGAGT ACCAGCACGC CTTCAAGACC CCGGATGCAG ATGCCGGTGA AGAAAGAGGA
AGCGGAGCTA CTAACCTCAG CCTGCTGAAG CAGGCTGGAG ACGTGGAGGA GAACCTGGA CCTATGACCG AGTACAAGCC
CACGGTGCGC CTCGCCACCC GCGACGAGT CCCCAGGGCC GTACGCACCC TCGCCGCGCG GTTCGCCGAC TACCCGCCA
CGCGCCACAC CGTCGATCCG GACCGCCACA TCGAGCGGT CACCGAGCTG CAAGAACTCT TCCTCACGCG CGTCGGGCTC
GACATCGGCA AGGTGTGGT CGCGGACGAC GCGCGCGCGG TGCGGTCTG GACCACGCCG GAGAGCGTCG AAGCGGGGGC
GGTGTTCGCC GAGATCGGCC CGCGCATGGC CGAGTTGAGC GGTTCGCCG TGCGCGCGCA GCAACAGATG GAAGGCCTCC
TGGCGCCGCA CCGGCCAAG GAGCCCGCT GGTTCCTGCG CACCGTCGGC GTCTCGCCCG ACCACAGGG CAAGGGTCTG
GGCAGCGCCG TCGTGCTCCC CGGAGTGGAG GCGGCCGAGC GCGCCGGGT GCCCGCCTTC CTGGAGACCT CCGCGCCCCG
CAACCTCCCC TTCTACGAGC GGCTCGGCTT CACCGTCACC GCGGACGTCG AGGTGCCGA AGGACCGCGC ACCTGGTGCA
TGACCCGCAA GCCCGGTGCC TGAAACTTGT TTATTGCAGC TTATAATGGT TACAAATAAA GCAATAGCAT CACAAATTC
ACAAATAAAG CATTTTTTTC ACTGCATTCT AGTTGTGGTT TGTCCAACT CATCAATGTA TCTTAATAAC TTCGTATAAT
GTATGCTATA CGAAGTTAT



OTI Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_001145468 , NM_001254732 , NM_001254733 , NM_015330
UniProt ID:	Q69YQ0
Synonyms:	CYTSA; GBBB2; OBLFC1
Summary:	This gene encodes a coiled-coil domain containing protein. The encoded protein may play a critical role in actin-cytoskeletal reorganization during facial morphogenesis. Mutations in this gene are a cause of oblique facial clefting-1. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. A read-through transcript composed of SPECCIL (sperm antigen with calponin homology and coiled-coil domains 1-like) and the downstream ADORA2A (adenosine A2a receptor) gene sequence has been identified, but it is thought to be non-coding. [provided by RefSeq, Jun 2013]

Product images:

