

Product datasheet for **KN502157**

Bhlhe22 Mouse Gene Knockout Kit (CRISPR)

Product data:

Product Type:	Knockout Kits (CRISPR)
Format:	2 gRNA vectors, 1 linear donor
Donor DNA:	EF1a-GFP-P2A-Puro
Symbol:	Bhlhe22
Locus ID:	59058



Components:

KN502157G1, Bhlhe22 gRNA vector 1 in pCas-Guide CRISPR vector (GEI00002)

KN502157G2, Bhlhe22 gRNA vector 2 in pCas-Guide CRISPR vector (GEI00002)

KN502157D, Linear donor DNA containing LoxP-EF1A-tGFP-P2A-Puro-LoxP:

The sequence below is cassette sequence only. The linear donor DNA also contains proprietary target sequence.

LoxP-EF1A-tGFP-P2A-Puro-LoxP (2739 bp)

ATAACTTCGT ATAATGTATG CTATACGAAG TTATCGTGAG GCTCCGGTGC CCGTCAGTGG GCAGAGCGCA CATCGCCAC
AGTCCCGCAG AAGTTGGGGG GAGGGGTCCG CAATTGAACC GGTGCCTAGA GAAGGTGGCG CGGGGTAAC TGGGAAAGTG
ATGTCGTGTA CTGGCTCCGC CTTTTTCCCG AGGGTGGGG AGAACCGTAT ATAAGTCAG TAGTCGCCGT GAACGTTCTT
TTTCGCAACG GGTTCGCCG CAGAACACAG GTAAGTGCCG TGTGTGGTTC CCGCGGGCCT GGCCTCTTTA CGGGTTATGG
CCCTTGCGTG CCTTGAATTA CTTCCACCTG GCTGCAGTAC GTGATTCTTG ATCCCGAGCT TCGGGTTGGA AGTGGGTGGG
AGAGTTCGAG GCCTTGCGCT TAAGGAGCCC CTTCGCCTCG TGCTTGAGTT GAGGCTGGC CTGGGCGCTG GGGCCGCCG
GTGCGAATCT GGTGGCACCT TCGCGCTGT CTCGCTGCT TCGATAAGTC TCTAGCCATT TAAATTTTTT GATGACCTGC
TGCGACGCTT TTTTCTGGC AAGATAGTCT TGTAATGCG GGCCAAAGAT TGCACACTGG TATTTCGGTT TTTGGGGCCG
CGGGCGGCGA CGGGGCCCGT GCGTCCCAGC GCACATGTTT GCGGAGGCGG GGCTGCGAG CGCGGCCACC GAGAATCGGA
CGGGGTAGT CTCAAGCTGG CCGGCCTGCT CTGGTGCCTG GCCTCGCGCC GCCGTGTATC GCGCCGCCCT GGGCGGCAAG
GCTGGCCCGG TCGGCACCAG TTGCGTGAGC GGAAGATGG CCGCTCCCG GCCCTGCTGC AGGGAGCTCA AAATGGAGGA
CGCGGCGCTC GGGAGAGCGG GCGGGTGAGT CACCCACACA AAGGAAAAGG GCCTTTCCTT CCTCAGCCGT CGCTTCATGT
GACTCCACGG AGTACCGGGC GCCGTCCAGG CACCTCGATT AGTTCCTGAG CTTTGGAGT ACCTCGTCTT TAGGTTGGGG
GGAGGGGTTT TATGCGATGG AGTTTCCCA CACTGAGTGG GTGGAGACTG AAGTTAGGCC AGCTTGGCAC TTGATGTAAT
TCTCCTGGA ATTTGCCCTT TTTGAGTTG GATCTTGTT CATTCTCAAG CTTGAGACAG TGGTTCAAAG TTTTCTCTT
CCATTTACAG TGTCGTGAAT GGAGAGCGAC GAGAGCGGCG TGCCCGCCAT GGAGATCGAG TGCCGCATCA CCGGCACCTC
GAACGCGCTG GAGTTCGAGC TGGTGGGCGG CGGAGAGGCG ACCCCGAGC AGGGCCGCAT GACCAACAAG ATGAAGAGCA
CCAAAGGCGC CCTGACCTTC AGCCCTACC TGCTGAGCCA CGTGATGGG TACGGCTTCT ACCACTTCGG CACTACCCC
AGCGGCTACG AGAACCCCTT CCTGCACGCC ATCAACAACG GCGGCTACAC CAACACCCG ATCGAGAAGT ACGAGGACGG
CGGCGTGCTG CACGTGAGCT TCAGCTACCG CTACGAGGCC GGCCGCGTGA TCGCGCACTT CAAGGTGATG GGCACCGGCT
TCCCCGAGGA CAGCGTGATC TTCACCGACA AGATCATCCG CAGCAACGCC ACCGTGGAGC ACCTGCACCC CATGGGCGAT
AACGATCTGG ATGGCAGCTT CACCCGCACC TTCAGCCTGC GCGACGGCGG CTAACAGC TCCGTGGTGG ACAGCCACAT
GCATTTCAAG AGCGCCATCC ACCCCAGCAT CCTGCAGAAC GGGGGCCCCA TGTTGCCTT CCGCCGCGTG GAGGAGGATC
ACAGCAACAC CGAGCTGGG ATCGTGAGT ACCAGCACGC CTTCAAGACC CCGGATGAG ATGCCGGTGA AGAAAGAGGA
AGCGGAGCTA CTAACCTCAG CCTGCTGAAG CAGGCTGGAG ACGTGGAGGA GAACCTGGA CCTATGACCG AGTACAAGCC
CACGGTGCGC CTCGCCACCC GCGACGAGT CCCCAGGGCC GTACGCACCC TCGCCGCGC GTTCGCCGAC TACCCGCCA
CGCGCCACAC CGTCGATCCG GACCGCCACA TCGAGCGGT CACCGAGCTG CAAGAACTCT TCCTCACGCG CGTCGGGCTC
GACATCGGCA AGGTGTGGT CGCGGACGAC GCGCGCGCGG TGCGGTCTG GACCACGCCG GAGAGCGTCG AAGCGGGGGC
GGTGTTCGCC GAGATCGGCC CGCGCATGGC CGAGTTGAGC GGTTCGCCG TGCCCGCGCA GCAACAGATG GAAGGCCTCC
TGGCGCCGCA CCGGCCAAG GAGCCCGCT GGTTCCTGCG CACCGTCGGC GTCTCGCCCG ACCACAGGG CAAGGGTCTG
GGCAGCGCCG TCGTGCTCCC CGGAGTGGAG GCGGCCGAGC GCGCCGGGT GCCCGCCTTC CTGGAGACCT CCGCGCCCCG
CAACCTCCCC TTCTACGAGC GGCTCGGCTT CACCGTCACC GCGGACGTCG AGGTGCCCGA AGGACCGCGC ACCTGGTGCA
TGACCCGCAA GCCCGGTGCC TGAAACTTGT TTATTGCAGC TTATAATGGT TACAAATAAA GCAATAGCAT CACAAATTC
ACAAATAAAG CATTTTTTTC ACTGCATTCT AGTTGTGGTT TGTCCAACT CATCAATGTA TCTTAATAAC TTCGTATAAT
GTATGCTATA CGAAGTTAT



OTI Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS. The kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process.
RefSeq:	NM_021560
UniProt ID:	Q8C6A8
Synonyms:	Beta3; Beta3a; Bhlhb5
Summary:	Inhibits DNA binding of TCF3/E47 homodimers and TCF3 (E47)/NEUROD1 heterodimers and acts as a strong repressor of Neurod1 and Myod-responsive genes, probably by heterodimerization with class a basic helix-loop-helix factors. Despite the presence of an intact basic domain, does not bind to DNA (By similarity). In the brain, may function as an area-specific transcription factor that regulates the postmitotic acquisition of area identities and elucidate the genetic hierarchy between progenitors and postmitotic neurons driving neocortical arealization. May be required for the survival of a specific population of inhibitory neurons in the superficial laminae of the spinal chord dorsal horn that may regulate pruritis. Seems to play a crucial role in the retinogenesis, in the specification of amacrine and bipolar subtypes. Forms with PRDM8 a transcriptional repressor complex controlling genes involved in neural development and neuronal differentiation (PubMed:22284184).[UniProtKB/Swiss-Prot Function]

Product images:

