


## Product datasheet for **LC403254**

### **FBXW7 (NM\_033632) Human Over-expression Lysate**


#### Product data:

Product Type:	Over-expression Lysates
Description:	FBXW7 HEK293T cell transient overexpression lysate (as WB positive control)
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC217398
Tag:	C-Myc/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<a href="#">NM_033632</a> , <a href="#">NP_361014</a>
Synonyms:	AGO; CDC4; FBW6; FBW7; FBX30; FBXO30; FBXW6; hAgo; hCdc4; SEL-10; SEL10
Predicted MW:	79.5 kDa
Components:	1 vial of 20 ug lyophilized gene specific transient over-expression cell lysate
Storage:	The lysate can be shipped at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C.
Preparation:	HEK293T cells in 10-cm dishes were transiently transfected with <a href="#">MegaTran</a> Transfection Reagent (TT200002) and 5ug <a href="#">TrueORF</a> cDNA plasmid. Transfected cells were cultured for 48hrs before collection. The cells were lysed in modified RIPA buffer (25mM Tris-HCl pH7.6, 150mM NaCl, 1% NP-40, 1mM EDTA, 1xProteinase inhibitor cocktail mix (Sigma), 1mM PMSF and 1mM Na3VO4), and then centrifuged to clarify the lysate. Protein concentration was measured by BCA kit (Thermo Scientific Inc.). To facilitate transportation and protein, the products are supplied as lyophilized proteins.
RefSeq:	<a href="#">NP_361014</a>
Locus ID:	55294
Cytogenetics:	4q31.3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Ubiquitin mediated proteolysis


[View online »](#)

**Product images:**


HEK293T cells were transfected with the pCMV6-Entry FBXW7 ([RC217398]) cDNA for 48 hrs and lysed. 5 ug of cell lysate was separated by SDS-PAGE and immunoblotted with anti-FBXW7 ([TA802883]). Positive lysates [LY403254] (100ug) and LC403254 (20ug) can be purchased separately from OriGene.



Western blot validation of overexpression lysate (Cat# [LY403254]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC217398] using transfection reagent MegaTran 2.0 (Cat# [TT210002]).