

Product datasheet for LC490001

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

SARS-CoV-2 N Protein (NC_045512) Virus Over-expression Lysate

Product data:

Product Type: Over-expression Lysates

Description: Transient overexpression lysate of SARS-CoV-2 nucleocapsid phosphoprotein (Severe acute

respiratory syndrome coronavirus 2)

Species: SARS-CoV-2

Expression cDNA Clone

or AA Sequence:

TrueORF Clone VC102555

Tag: C-MYC/DDK

Detection Antibodies: Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)

ACCN: NC 045512, YP 009724397

Predicted MW: 45.6 kDa

Components: 1 vial of 20 ug lyophilized gene specific transient over-expression cell lysate

Storage: The lysate can be shipped at ambient temperature. Upon receiving, store the sample at -

20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt

when stored at -20°C.

Bioactivity: N/A

Preparation: HEK293T cells in 10-cm dishes were transiently transfected with MegaTran Transfection

Reagent (TT200002) and 5ug TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs before collection. The cells were lysed in modified RIPA buffer (25mM Tris-HCl pH7.6, 150mM NaCl, 1% NP-40, 1mM EDTA, 1xProteinase inhibitor cocktail mix (Sigma), 1mM PMSF and 1mM Na3VO4), and then centrifuged to clarify the lysate. Protein concentration was measured by BCA kit (Thermo Scientific Inc.). To facilitate transportation and protein, the


products are supplied as lyophilized proteins.

RefSeq: <u>YP 009724397</u>

Locus ID: 43740578



Product images:



Western blot validation of overexpression lysate (Cat# [LY490001]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [VC102555] using transfection reagent MegaTran 2.0 (Cat# [TT210002])