

## Product datasheet for **LC705973**

### **B4galt6 (NM\_019737) Mouse Over-expression Lysate**


#### **Product data:**

<b>Product Type:</b>	Over-expression Lysates
<b>Description:</b>	B4galt6 HEK293T cell transient overexpression lysate (as WB positive control)
<b>Species:</b>	Mouse
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone MR205973
<b>Tag:</b>	C-Myc/DDK
<b>Detection Antibodies:</b>	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
<b>ACCN:</b>	<u><a href="#">NM_019737</a></u> , <u><a href="#">NP_062711</a></u>
<b>Synonyms:</b>	AA536803; AU022389
<b>Predicted MW:</b>	44.8 kDa
<b>Components:</b>	1 vial of 20 ug lyophilized gene specific transient over-expression cell lysate
<b>Storage:</b>	The lysate can be shipped at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples are stable for 12 months from date of receipt when stored at -20°C. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples can be reconstituted with SDS Sample Buffer. Prior to SDS-PAGE fractionation, boil the lysate for 5 minutes. Lysate samples are stable for 12 months from date of receipt when stored at -20°C.
<b>Preparation:</b>	HEK293T cells in 10-cm dishes were transiently transfected with <u><a href="#">MegaTran</a></u> Transfection Reagent (TT200002) and 5ug <u><a href="#">TrueORE</a></u> cDNA plasmid. Transfected cells were cultured for 48hrs before collection. The cells were lysed in modified RIPA buffer (25mM Tris-HCl pH7.6, 150mM NaCl, 1% NP-40, 1mM EDTA, 1xProteinase inhibitor cocktail mix (Sigma), 1mM PMSF and 1mM Na3VO4), and then centrifuged to clarify the lysate. Protein concentration was measured by BCA kit (Thermo Scientific Inc.). To facilitate transportation and protein, the products are supplied as lyophilized proteins.
<b>RefSeq:</b>	<u><a href="#">NP_062711</a></u>
<b>Locus ID:</b>	56386
<b>Cytogenetics:</b>	18 A2



[View online »](#)

## Product images:



Western validation with an anti-DDK antibody; L: Control HEK293 lysate R: Over-expression lysate