

Product datasheet for LY411500

FOXL2 (NM_023067) Human Over-expression Lysate

Product data:

Product Type: Over-expression Lysates Description: Transient overexpression lysate of forkhead box L2 (FOXL2) Species: Human HEK293T **Expression Host:** TrueORF Clone RC209127 **Expression cDNA Clone** or AA Sequence: Tag: C-Myc/DDK **Detection Antibodies:** Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100) ACCN: NM 023067, NP 075555 BPES; BPES1; PFRK; PINTO; POF3 Synonyms: Predicted MW: 38.6 kDa 1 vial of 100 µg gene specific transient over-expression cell lysate in RIPA buffer **Components:** 1 vial of 100 µg whole HEK293T cell lysate in RIPA buffer 1 vial of 250ul 2xSDS Sample Buffer (4% SDS, 125mM Tris-HCl pH6.8, 10% Glycerol, 0.002% Bromophenol blue, 100mM DTT) The lysate is shipped with dry ice. Upon receiving, store the sample at -80°C. Also after Storage: dilution, the protein sample should be aliquoted and stored at -80°C for long term storage. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C. **Preparation:** HEK293T cells in 10-cm dishes were transiently transfected with MegaTran Transfection Reagent (TT200002) and 5ug TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs before collection. The cells were lysed in modified RIPA buffer (25mM Tris-HCl pH7.6, 150mM NaCl, 1% NP-40, 1mM EDTA, 1xProteinase inhibitor cocktail mix (Sigma), 1mM PMSF and 1mM Na3VO4), and then centrifuged to clarify the lysate. Protein concentration was measured by BCA kit (Thermo Scientific Inc.). Cell lysates were aliquoted and stored at -20°C before shipping. RefSeq: NP 075555 Locus ID: 668



/iew online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Sourigene FOXL2 (NM_023067) Human Over-expression Lysate – LY411500


Cytogenetics:

3q22.3

Protein Families:

Druggable Genome, Transcription Factors

Product images:



Western blot validation of overexpression lysate (Cat# LY411500) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC209127] using transfection reagent MegaTran 2.0 (Cat# [TT210002]).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US