

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for LY421213

#### ZNF133 (NM\_001083330) Human Over-expression Lysate

#### **Product data:**

Description:Transient overexpression lysate of zinc finger protein 133 (ZNF133), transcript variant 2Species:HumanExpression Host:HEK293TFxpression cDNA CloneTrueORF Clone RC223016	
Expression Host: HEK293T	
Expression cDNA Clone TrueORF Clone RC223016	
or AA Sequence:	
Tag: C-Myc/DDK	
<b>Detection Antibodies:</b> Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)	
ACCN: <u>NM 001083330</u> , <u>NP 001076799</u>	
Synonyms: pHZ-13; pHZ-66; ZNF150	
Predicted MW: 73.3 kDa	
Components:1 vial of 100 μg gene specific transient over-expression cell lysate in RIPA buffer1 vial of 100 μg whole HEK293T cell lysate in RIPA buffer1 vial of 250ul 2xSDS Sample Buffer (4% SDS, 125mM Tris-HCl pH6.8, 10% Glycerol, 0.002%Bromophenol blue, 100mM DTT)	
Storage:The lysate is shipped with dry ice. Upon receiving, store the sample at -80°C. Also after dilution, the protein sample should be aliquoted and stored at -80°C for long term storage Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffe provided. Lysate samples are stable for 12 months from the date of receipt when stored at 80°C.	er
Preparation:HEK293T cells in 10-cm dishes were transiently transfected with MegaTran Transfection Reagent (TT200002) and 5ug TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs before collection. The cells were lysed in modified RIPA buffer (25mM Tris-HCl pH7.6 150mM NaCl, 1% NP-40, 1mM EDTA, 1xProteinase inhibitor cocktail mix (Sigma), 1mM PMS and 1mM Na3VO4), and then centrifuged to clarify the lysate. Protein concentration was measured by BCA kit (Thermo Scientific Inc.). Cell lysates were aliquoted and stored at -200 before shipping.	F
RefSeq: NP 001076799	
Locus ID: 7692	



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US


## STATE STATE

Cytogenetics:

20p11.23

Protein Families: Transcription Factors

### **Product images:**



Western blot validation of overexpression lysate (Cat# LY421213) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC223016] using transfection reagent MegaTran 2.0 (Cat# [TT210002]).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US