

Product datasheet for **LY438876**


GLT25D2 (COLGALT2) (NM_001303420) Human Over-expression Lysate

Product data:

Product Type:	Over-expression Lysates
Description:	Transient overexpression lysate of collagen beta(1-O)galactosyltransferase 2 (COLGALT2), transcript variant 2
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238876
Tag:	C-MYC/DDK
ACCN:	<u>NM_001303420</u> , <u>NP_001290349</u>
Synonyms:	C1orf17; ColGalT 2; GLT25D2
Predicted MW:	65.1 kDa
Components:	1 vial of 100 µg gene specific transient over-expression cell lysate in RIPA buffer 1 vial of 100 µg whole HEK293T cell lysate in RIPA buffer 1 vial of 250ul 2xSDS Sample Buffer (4% SDS, 125mM Tris-HCl pH6.8, 10% Glycerol, 0.002% Bromophenol blue, 100mM DTT)
Storage:	The lysate is shipped with dry ice. Upon receiving, store the sample at -80°C. Avoid repeated freeze-thaw cycles. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C.
Preparation:	HEK293T cells in 10-cm dishes were transiently transfected with PEI and 5ug TrueORF cDNA plasmid. Transfected cells were cultured for 72hrs before collection. The cells were lysed in modified RIPA buffer (25mM Tris-HCl pH7.6, 150mM NaCl, 1% NP-40, 1mM EDTA, 1xProteinase inhibitor cocktail mix (Sigma), 1mM PMSF and 1mM Na3VO4), and then centrifuged to clarify the lysate. Protein concentration was measured by BCA kit (Thermo Scientific Inc.).
RefSeq:	<u>NP_001290349</u>
Locus ID:	23127
Cytogenetics:	1q25.3
Protein Families:	WB


[View online »](#)

Product images:



Western blot validation of overexpression lysate (Cat# LY438876) using anti-DDK antibody (Cat# [TA592569]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC238876] using transfection reagent PEI.