

## Product datasheet for **LY500031**

### IGROV1 (Human Ovarian Carcinoma) Whole Cell Lysate

#### Product data:

Product Type:	Cell Line Lysates
Tissue:	Ovary
Host:	Human
Lysis Buffer:	Modified RIPA buffer (25mM Tris-HCl pH7.6, 150mM NaCl, 1% NP-40, 1mM EDTA, 1xProteinase inhibitor cocktail mix (Sigma) and 1mM Na <sub>3</sub> VO <sub>4</sub> )
Components:	1 vial of 100ug whole cell lysate in RIPA buffer 1 vial of 250ul 2xSDS Sample Buffer (4% SDS, 125mM Tris-HCl pH6.8, 10% Glycerol, 0.002% Bromphenol blue, 100mM DTT)
Recommended Applications:	IP, WB
Storage:	Shipped on dry ice. Store the lysate at -80°C and avoid repeated freeze / thaw cycles. The lysate is stable for at least one year from date of shipping when stored at -80°C. After incubation with SDS-PAGE sample buffer, SDS-PAGE samples need to be stored in -20°C for long term storage.



[View online »](#)