

Product datasheet for **LY500065**

SV-T2 (Mouse BALB/3T3 embryo fibroblast cell, SV40 transformed, clone A31) Whole Cell Lysate

Product data:

Product Type:	Cell Line Lysates
Tissue:	Embryo
Host:	Mouse
Lysis Buffer:	Modified RIPA buffer (25mM Tris-HCl pH7.6, 150mM NaCl, 1% NP-40, 1mM EDTA, 1xProteinase inhibitor cocktail mix (Sigma) and 1mM Na ₃ VO ₄)
Components:	1 vial of 100ug whole cell lysate in RIPA buffer 1 vial of 250ul 2xSDS Sample Buffer (4% SDS, 125mM Tris-HCl pH6.8, 10% Glycerol, 0.002% Bromphenol blue, 100mM DTT)
Recommended Applications:	IP, WB
Storage:	Shipped on dry ice. Store the lysate at -80°C and avoid repeated freeze / thaw cycles. The lysate is stable for at least one year from date of shipping when stored at -80°C. After incubation with SDS-PAGE sample buffer, SDS-PAGE samples need to be stored in -20°C for long term storage.



[View online »](#)