

Product datasheet for **MC200003**

Snrpa (NM_001046637) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Snrpa (NM_001046637) Mouse Untagged Clone
Tag: Tag Free
Symbol: Snrpa
Synonyms: C430021M15Rik; Rnu1a-1; Rnu1a1; U1-A; U1A
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC003229 sequence for NM_001046637
 GGGGGCGGGACCCCAAGTGTTCCTTTGGCAGACGCTGAGGCCGGGGGAGAGACTCTGAACCTACCGAGTT
 TGCCCGCGTGATTAGCTCGCGGAAGGGAGAAGCCGTTGCCGGGAAGATAGTTGTGTGTCTGTGCAGCG
 TAGCTCGGCGGGCGTTGTGCGCGGGAGCGACGGCTCAGGTGCCGTACCTGAACATCTCACACCATGGCCA
 CCATAGCCACCATGCCAGTTCGCCGAGACCCGTGCCAACCACTATTTATATCAACAATCTCAATGAGAA
 GATCAAGAAGGATGAGCTCAAGAAGTCCCTGTATGCCATCTTCTCCAGTTTGGCCAGATCCTGGATATC
 CTGGTGTCTCGGATCATGAAGATGAGGGGCCAGGCCTTCGTATCTTCAAGGAGGTCAACAGCGCCACCA
 ATGCCCTGCGCTCCATGCAGGGCTTCCCTTTCTACGACAAGCCCATGCGCATCCAGTACGCAAAGACTGA
 CTCGGACATCATTGCCAAGATGAAGGGCACCTATGTGGAGAGAGATCGCAAACGAGAGAAGAGGAAGCCC
 AAGAGTCAGGAGACACCTGCTGCCAAAAAGGCTGTTCAAGGTGGGGCAGCTGCGCCCGTGGTGGGGGCTG
 TCCAGCCCGTACCGGGCATGCCACCGATGCCTCAGGCACCCCGCATCATGCACCATATGCCAGGACAGCC
 TCCTACATGCCGCCACCTGGCATGATCCCGCCACCAGGCCTCGCTCCTGGCCAGATCCCCCTGGGGCC
 ATGCCCCCACAGCAGCTCATGCCTGGGCAGATGCCGCCTGCCAGCCTCTCTCCGAGAACCACCCAATC
 ACATCCTGTTCTCACCAACCTGCCTGAGGAGACCAACGAGCTCATGCTCTCCATGCTCTTCAACCAAGTT
 CCCTGGCTTCAAGGAGGTGCGTCTGGTCCCTGGGGCCCATGACATCGCCTTCGTGGAGTTTGACAATGAA
 GTGCAGGCTGGGGCAGCACGAGATGCCCTGCAAGGCTTTAAGATCACACAAAACAATGCTATGAAGATCT
 CTTTTGCCAAGAAGTAGCGCCTTCCCTATGGAGTACCCAGTCCCTTCCCCCTCCCTTGCTCAGTC
 CCTGAAGGTAAGTCCCCCTTAGGGACCTTCTCAGAGCCGTGTGTGTGGTTGCCACACAGCATTGTACC
 CAGAGTCTGTCCCAGACATTGCACCTGGCGCTGTAGATTGTGATTAAGTGAGTTTTTGAGGTTTGGT
 TTTTTAAAAAAAAAAAAAAAAA

Restriction Sites: RsrII-NotI
ACCN: NM_001046637
Insert Size: 864 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC003229 , AAH03229
RefSeq Size:	1280 bp
RefSeq ORF:	864 bp
Locus ID:	53607
UniProt ID:	Q62189
Cytogenetics:	7 A3
Gene Summary:	<p>Component of the spliceosomal U1 snRNP, which is essential for recognition of the pre-mRNA 5' splice-site and the subsequent assembly of the spliceosome. U1 snRNP is the first snRNP to interact with pre-mRNA. This interaction is required for the subsequent binding of U2 snRNP and the U4/U6/U5 tri-snRNP. SNRPA binds stem loop II of U1 snRNA. In a snRNP-free form (SF-A) may be involved in coupled pre-mRNA splicing and polyadenylation process. May bind preferentially to the 5'-UGCAC-3' motif on RNAs (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 3 encode the same protein.</p>